

Tutorial:Introduction to Cytoscape

Slideshow Introduction to Cytoscape (about 30 minutes) ^[1]

Handout Introduction_to_Cytoscape.pdf (16 pages) ^[2]

Tutorial Sources Cytoscape Tutorial (Yeyejide Adeleye) ^[3]

Tutorial Curators Anna Kuchinsky, Scooter Morris, Alex Pico

Cytoscape is an open source software tool for *integrating*, *visualizing*, and *analyzing* data in the context of networks. This tutorial will cover:

1. Navigating Cytoscape
2. Visualizing Data on Networks


Part two of this tutorial can be found [here](#)

Navigating Cytoscape

This section will introduce the Cytoscape user interface. First of all we will look at the basic UI of Cytoscape. Then we will show all menu features of Cytoscape and the extended functionality provided by plugins.

Cytoscape Layout and User Interface

Launch Cytoscape. You should see a window that looks like this:



- At the top of the Cytoscape Desktop window is the toolbar, which contains the command buttons. The name of each command button is shown when the mouse pointer hovers over it.
- In the upper right is the Main Network View window, where network data will be displayed. This region is initially blank.
- At left is the Control Panel (Network Management) Panel. This lists the available networks by name and provides information on the number of nodes and edges.

- Immediately below the Control Panel is the Network Overview Pane
- At lower right is the Data Panel which can be used to display node, edge, and network attribute data

The Network Management and Data browser panels are dockable tabbed panels known as CytoPanels. You can undock any of these panels by clicking on the Float Window control in the upper-right corner of the CytoPanel. The Data Panel starts off with three tabs: **Node Attribute Browser**, **Edge Attribute Browser**, and **Network Attribute Browser**; the Network Management panel starts off with four tabs: **Network**, **VizMapper**, **Editor**, and **Filters**. Loaded plugins might add tabs to either of these CytoPanels.

Cytoscape Menus

We will briefly run through all the menus available in Cytoscape.

File

The File menu contains basic file functionality:

- **File → Open** for opening a Cytoscape session file
- **File → New** for creating a new network
- **File → Save** for saving a session file
- **File → Import** for importing data such as networks and attributes
- **File → Export** for exporting data and images.
- **File → Print** allows printing
- **File → Quit** closes all windows of Cytoscape and exits the program

Edit

The Edit menu contains:

- Undo and Redo functions which undo and redo edits made in the Attribute Browser, the Network Editor and the Layout.
- Options for creating and destroying views (graphical representations of a network) and networks
- Options for deleting selected nodes and edges from the current network.
- All deleted nodes and edges can be restored to the network via **Edit → Undo**.
- **Edit → Preferences → Properties** to edit preferences for properties and plugins

View

The View menu allows you to display or hide:

- The network management panel (Control Panel)
- The attribute browser (Data Panel)
- Results Panel
- VizMapper

Select

The Select menu contains:

- Options for selecting nodes and edges
- The **Select → Use Filters** option allows filters to be created for automatic selection of portions of a network whose node or edge attributes meet a filtering criterion (see below for the filters section).

Layout

The Layout menu has an array of features for visually organizing the network:

- Rotate, Scale, Align and Distribute are tools for manipulating the network visualization.
- The bottom section of the menu lists a variety of layout algorithms which automatically lay a network out.

Plugins

The Plugins menu contains options for managing your plugins (install/update/delete) and may have options added by plugins that have already been installed, such as the Agilent Literature Search or Merge Networks.

- Depending on which plugins are loaded, the plugins that you see may be different than what appear here.


Help

- The Help menu allows you to launch the help viewer and browse the table of contents for this manual.
- The “About...” option displays information about the running version of Cytoscape.

Loading a Simple Network

- Go to **File-> Import -> Network (multiple file types)**
- You should see the Import Network File Dialog
- For Data Source Type select **Local** and then click **Select**
- Open the **sampleData** folder and select **galFiltered.sif** and then click on **Open** and then **Import**



You should see the following:



The SIF file format is about as simple as it gets. It consists of 3 columns: source, interaction type, and target. “Source” and “target” are gene/protein identifiers that are used to define nodes, while “interaction type” serves to label the edge connecting each pair of nodes.

Manipulating Your Network

Now that you have a network loaded, you can interact with it in a number of ways:

- Start by clicking on the node at the upper left corner of the network. The node will turn yellow. If you hold your mouse down over the node and drag it around the node will move on the screen.
- Now add another node to the selection by holding down the Shift key and clicking on a node. Note that both nodes are now selected (yellow). Again, move the nodes around. Note that both nodes will move.
- To select a group of nodes, hold the mouse down in the upper left-hand corner and drag your mouse over a region of the network. Again, a group of nodes will be selected and can be moved around on the screen.
- To zoom in on the selected nodes, click on the  icon.
- To move the window around the network, you can either use the middle mouse button, or drag the small window outlined in blue around in the Network Overview Pane.
- Finally, zoom your network out by clicking on the  icon.

While useful, hand selecting nodes in dense networks can be error-prone and difficult. However, you can specifically search for a node by name or attribute:

- In the **Search:** box at the top of the screen, type in *ynr050c*. This will select that node and zoom the display to focus on it.

The **Search:** box will also allow you to select nodes by other attributes, but first, we need to import more attributes...

Visualizing Data on Networks


Cytoscape provides a number of features to load arbitrary data and visualize that data by mapping the attribute values to visual styles.

Importing Your Data

Cytoscape can read file structures that are delimited Text or Excel files.

- Go to **File-> Import->Attribute from Table (Text/MS Excel)**
- You should see the **Import Attribute from Table** Dialog Box
- Make sure the **Node** radio attribute button is selected in **Data Sources Attributes**
- In **Data Sources**, click on **Select File**
- Select **galExpData.pvals** from the **sampleData** folder and then click on **Open**
- In the **Advanced** section, select **Show Text File Import Options**
- By default, **Tab** is selected in the **Delimiter** section. Instead, select **Space**.
- In the **Attribute Names** section, select **Transfer first line attribute names**
- Also in the **Advanced** section, select **Import Everything**


*The 'Import Everything' option tells Cytoscape to load **all** of your data, not just the records that match currently loaded networks.*



If you were to click 'Import' now you would see a pop-up message complaining about duplicate attribute names. Notice that columns 6, 7 and 8 have the same names as 3, 4 and 5. The next steps show you how to fix column names without having to start over from Excel.

- Right click on the column header of the first duplicate name (column 6, gal1RG) to open the **Set Attribute Name and Type** dialog box
- Add the suffix 'sig' to the name (e.g., 'gal1RGsig') to distinguish the column as containing p-values (significance)
- Repeat these last two steps on columns 7 and 8 (gal4RG and gal80R)
- Now click **Import**

You've successfully loaded data! Now explore the **Data Panel** to confirm the mapping of the data to the network.

- Locate the **Select Attributes** button  in the **Data Panel**
- Select the data attributes: gal1RG, gal4RG, and gal80R
- Return to the **Data Panel** by click away from the attribute list
- Select some nodes in the network (ctrl-A selects all nodes) and see the associated data in the **Data Panel**

You now have access to your networked data and can begin playing with visualization.


- Go to **Layout -> Cytoscape Layouts -> Force-Directed Layout**
- Go to the **VizMapper** tab in the **Control Panel** and select **Sample 1**


These are default visual styles. In the next section you will customize a visual style to highlight your data values.

Viewing Attributes

In the previous section, you were able to select a couple of attributes to display in the Data Panel. In this section, we will explore a little more about attributes and the Data Panel.

- First click on the **Configure search options** icon: .
- Change the **Select Attribute:** to **COMMON**. You should see this screen:



- Click on **Apply**
- Now type *mcm1* in the **Search:** box. This will select the node: *YMR043W*, and display the attributes for that node in the **Data Panel**.
- We're going to add a new attribute for *MCM1*. Click on the  icon and select **String Attribute**.
- Type in **pdb** for the name of the attribute -- this will define a new string attribute for nodes, and add it to the **Data Panel**.
- Now click into the empty cell for newly-created **pdb** attribute for **YMR043W** and type *1mm*, which is the PDB ID for the yeast protein *mcm1*. You need to hit Return or Tab to enter that data.
- Move the **pdb** attribute to be the second column by dragging the column header to be behind the **ID** column
- Finally, select a number of nodes and note that the attributes for all of the nodes are shown in the **Data Panel**
- By clicking on the column header, you can sort the columns. Clicking again changes the order of the sort.

Visualizing Data with VizMapper


- Go to **File-> Open** (click 'Yes' to losing current session)
- You should see the Open a Session File Dialog
- Open the **sampleData** folder and select **galFiltered.cys** and then click on **Open**

Notice how the galFiltered Style maps multiple data and annotation values to:


- Edge Color
- Edge Line Style
- Node Border
- Node Color
- Node Label
- Node Size
- Node Tooltip

Modifying a Visual Style

Customising the way you visualize and manipulate networks is a key function of Cytoscape. This is achieved through the use of the VizMapper tool.

- To launch the VizMapper, either select the VizMapper tab on the Control Panel or click on the VizMapper icon  at the top of the tool bar
- Find **Node Color** in the **Visual Mapping Browser** and expand it by clicking on the triangle icon for expand/collapse
- Click on the value 'gal4RGexp' to select an alternative data value to map to node color: select 'gal80Rexp'. Notice the changes in the network display.
- Click on the gradient color mapping to open the **Gradient Editor**
- Control the values and colors of the mapping by means of triangular handles and endpoint markers
- Double-click on any handle or marker to change it's color: change green to blue and change red to yellow. Notice the immediate changes to the network
- Click-and-drag any handle to slide its value between the min and max
- To save your changes, close the editor
- Find **Edge Color** in the **Visual Mapping Browser** and expand it by clicking on the triangle icon
- Notice the discrete mapping of colors to values. Click on any color mapping to change the color.
- Close the editor to save your changes


Creating a Visual Style (Advanced)


- In the **Current Visual Style** section, start a new visual style by clicking the **Options** button  and selecting "Create new Visual Style"
- Enter a name for your custom visual style when prompted, click on **OK**
- Find **Node Color** in the list of **Unused Properties** and double-click to activate
- For **Node Color** select the value 'gal1RGexp' to map expression fold values. For **Mapping Type** select 'Continuous'. For **Graphical View** click on the gradient to open the dialog:
 - Click on **Min/Max** button and set to -1 and 1, respectively
 - Double-click on gradient handles to set colors, e.g.
 - below -1.0 = green
 - -0.8 = green
 - 0.8 = red
 - above 1.0 = red
 - Click **Add** to add another gradient handle (added to max by default). Drag to center at 0.0. Leave as white color.
 - Close gradient dialog
- Next, activate **Node Size**. Select 'gal1RGsig' to map p-values. For **Mapping Type** select 'Continuous'. Click on gradient to open dialog:
 - Note: We want smaller p-values (more significant) to show as larger nodes
 - Click **Min/Max** to verify the correct values: 1.0E-8 and 1.0E-3, respectively.
 - Double-click on solid red box for p-values *below the minimum* (far left) to set the *maximum* node size. Set to 70.0
 - Double-click on solid red box for p-values *below the maximum* (far right) to set the *minimum* node size. Set to 20.0
 - Select handles along the top border to drag the x-position (p-value) of the gradient points. When selected, you can also set the Attribute and Node Size value by editing the **Handle Settings** field. You can also drag (or

double-click) the open red boxes to set the y-position (Node Size).

- Set the left open box to Node Size 70.0 and the right open box to 20.0. This defines the range of the gradient.
- Click 'Add' to add a new gradient handle. Set it to Attribute Value = 1E-4 and Node Size = 40.0.

This creates a pseudo-exponential gradient mapping:



- Click 'OK' in the gradient dialog and **explore** the visualization you've created!
- Zoom by selecting an area and clicking 
- Use the bird's-eye-view panel (bottom-left) to pan around network
- Return to VizMapper and switch mappings to another column of data:
 - For **Node Color** select 'gal4RG'
 - For **Node Size** select 'gal4RGsig'

Notice the change in the view of the network! You can reuse the mappings across multiple datasets. Now is a good time to save your Cytoscape session, to save the visual style you've created.

Laying Out Your Network


A network layout is a process that positions the nodes and edges for the network. There are a large variety of layouts in Cytoscape and plugins might add new layouts. All of the layouts will appear under the **Layouts** menu. In this section, we'll explore some of the layouts in the **Cytoscape Layouts** category, which are the core layouts supported by the Cytoscape team. All of the **Cytoscape Layouts** support the ability to only layout a portion of the network, and most expose parameters that can be used to tune the layout algorithm.

Simple Layouts

Grid Layout

The simplest layout that Cytoscape provides is the **Grid Layout**, which simply places all of the nodes in a grid arrangement.

- Using the network you loaded before, select **Layout->Cytoscape Layouts->Grid**. You should see the image below:




- Now, select **Edit->Undo Grid Layout Layout** and notice how your network reverts to the previous condition.
 - Many (but not all) of Cytoscape's actions may be undone.
- Select a group of nodes in your network and again select **Layout->Cytoscape Layouts->Grid**.
 - Notice that this now has a sub-menu with two options: **All Nodes** and **Selected Nodes Only**.
- Choose **Selected Nodes Only**.
 - Note that only the nodes you selected are changed.

Circular Layout


Grid is very fast, but often not very helpful. Alternatives:

- Circular Layout** places all of the nodes in a circular arrangement.
 - Very quick
 - Usually not very informative.
 - Partitions* the network into disconnected parts and independently lays out those parts.



Hierarchical Layout

- **Hierarchical Layout** forces the nodes into a tree structure.
 - Works best when the network is naturally tree-structured
 - Also works reasonably well when the network is mostly hierarchical.



Data-Driven Simple Layouts


Often what is desired is to organize the nodes in space to reflect some data property of the nodes themselves. We'll look at three of the simple ones:

- **Degree Sorted Circle:** Orders the node around a circle based on node degree (number of edges)
- **Attribute Circule Layout:** Orders the node around a circle based on the value of some attribute
- **Group Attributes Layout:** Groups the nodes based on the value of some attribute

Degree Sorted Circle Layout


For example, assume you are interested in hubs (nodes with high degree):

- Select **Layout->Degree Sorted Circle Layout**.
 - Note the highest degree nodes are in the same region of the circle and the degree decreases as you proceed counter-clockwise around the circle.
 - Also note that this layout supports partitioning of the network into disconnected components.




Attribute Circle Layout

- Now try **Layout->Attribute Circle Layout->Degree**.
 - This should give you a very similar layout to the Degree sorted circle layout.
 - You should also try using other attributes to layout.



Group Attributes Layout

- Finally, select **Layout->Cytoscape Layouts->Group Attributes Layout->Degree**.
 - Note that this layout organizes all of the nodes with the same degree into a circle and then positions the circles into a grid.
 - Partitions the graph before layouts.



These layouts are useful to represent data attributes associated with the nodes of a network, however, they don't provide much information about the network structure itself. For that, using more complicated layouts

are required.

Force-Directed Layouts



The force-directed methods for laying out graphs all use some kind of physical simulation that models the nodes as physical objects and the edges as springs connecting those objects together.

- Cytoscape provides 4 force-directed layouts:
 - Edge-weighted Force directed (BioLayout)
 - Edge-weighted Spring Embedded
 - Spring Embedded
 - Force-Directed

In this exercise, we'll use the **Force-Directed Layout**, which is a port of the layout by the same name in the *prefuse* package (see [<http://prefuse.org>]^[4]).


Force-Directed Layout

- We'll do all of our work from the Layout settings dialog, so bring up the dialog by **Layout->Settings...**
- In the **Layout Settings** dialog:
 - Select **Force-Directed Layout** under **Select algorithm to view settings**.
 - The result should look like the image on the right
- Lets start by creating a layout using the default parameters.
 - Select **Execute Layout**.
 - Note how this exposes more of the structure of the network than any of the previous layouts.




Now, let's see what the parameters do.

- Start by changing the **Default Spring Length** to **20**, and select **Execute Layout**.
- The default spring length is the length of a spring (edge) with no forces exerted upon it. Essentially, it's the length of an edge connecting two nodes with no other connections.
- Note how the layout has changed.
- The resulting network tends to be more closely packed in the denser regions.



- Finally, let's change **Default Spring Length** back to **50** and the **Default Spring Coefficient** to **1E-3**.
- The spring coefficient is a measure of the strength of the spring.
- Now click **Execute Layout**.
- Note that the network remains pretty compact, even though we've reset the default spring length back to 50.



Force-Directed Layout is also a good layout to use for circumstances when you want to have the length of your edges reflect some numeric weight on each edge.

- Select the attribute containing your edge weights from the menu **The edge attribute that contains the weights**.
- This often requires some tuning of the spring coefficient and the spring length to get an aesthetically pleasing layout.

References

- [1] http://opentutorials.rbvi.ucsf.edu/index.php?title=Tutorial:Introduction_to_Cytoscape&ce_slide=true&ce_style=cytoscape
 - [2] http://opentutorials.rbvi.ucsf.edu/index.php/File:Introduction_to_Cytoscape.pdf
 - [3] <http://cytoscape.wodaklab.org/wiki/CytoRetreat2007Tutorials>
 - [4] <http://prefuse.org>
-

Article Sources and Contributors

Tutorial:Introduction to Cytoscape *Source:* <http://opentutorials.cgl.ucsf.edu/index.php?oldid=2161> *Contributors:* AlexanderPico, AnnaKuchinsky, KristinaHanspers, ScooterMorris

Image Sources, Licenses and Contributors

Image:screen.png *Source:* <http://opentutorials.cgl.ucsf.edu/index.php?title=File:Screen.png> *License:* unknown *Contributors:* AlexanderPico, AnnaKuchinsky

Image:loadingnetwork.png *Source:* <http://opentutorials.cgl.ucsf.edu/index.php?title=File>Loadingnetwork.png> *License:* unknown *Contributors:* AlexanderPico, AnnaKuchinsky

File:Zoom.png *Source:* <http://opentutorials.cgl.ucsf.edu/index.php?title=File:Zoom.png> *License:* unknown *Contributors:* ScooterMorris

File:ZoomOut.png *Source:* <http://opentutorials.cgl.ucsf.edu/index.php?title=File:ZoomOut.png> *License:* unknown *Contributors:* ScooterMorris

Image:importing.png *Source:* <http://opentutorials.cgl.ucsf.edu/index.php?title=File:Importing.png> *License:* unknown *Contributors:* AnnaKuchinsky

Image:attributebutton.png *Source:* <http://opentutorials.cgl.ucsf.edu/index.php?title=File:Attributebutton.png> *License:* unknown *Contributors:* AnnaKuchinsky

File:QuickFind.png *Source:* <http://opentutorials.cgl.ucsf.edu/index.php?title=File:QuickFind.png> *License:* unknown *Contributors:* ScooterMorris

Image:QuickFindConfiguration.png *Source:* <http://opentutorials.cgl.ucsf.edu/index.php?title=File:QuickFindConfiguration.png> *License:* unknown *Contributors:* ScooterMorris

File:NewAttribute.png *Source:* <http://opentutorials.cgl.ucsf.edu/index.php?title=File:NewAttribute.png> *License:* unknown *Contributors:* ScooterMorris

Image:vizmapper.png *Source:* <http://opentutorials.cgl.ucsf.edu/index.php?title=File:Vizmapper.png> *License:* unknown *Contributors:* AlexanderPico, AnnaKuchinsky, KristinaHanspers

Image:visualstyle.png *Source:* <http://opentutorials.cgl.ucsf.edu/index.php?title=File:Visualstyle.png> *License:* unknown *Contributors:* AnnaKuchinsky

Image:gradient.png *Source:* <http://opentutorials.cgl.ucsf.edu/index.php?title=File:Gradient.png> *License:* unknown *Contributors:* AlexanderPico, ScooterMorris

Image:zoombutton.png *Source:* <http://opentutorials.cgl.ucsf.edu/index.php?title=File:Zoombutton.png> *License:* unknown *Contributors:* AnnaKuchinsky

Image:GridLayout.png *Source:* <http://opentutorials.cgl.ucsf.edu/index.php?title=File:GridLayout.png> *License:* unknown *Contributors:* ScooterMorris

Image:CircularLayout.png *Source:* <http://opentutorials.cgl.ucsf.edu/index.php?title=File:CircularLayout.png> *License:* unknown *Contributors:* ScooterMorris

Image:HierarchicalLayout.png *Source:* <http://opentutorials.cgl.ucsf.edu/index.php?title=File:HierarchicalLayout.png> *License:* unknown *Contributors:* ScooterMorris

Image:DegreeSortedLayout.png *Source:* <http://opentutorials.cgl.ucsf.edu/index.php?title=File:DegreeSortedLayout.png> *License:* unknown *Contributors:* ScooterMorris

Image:GroupAttributesLayout.png *Source:* <http://opentutorials.cgl.ucsf.edu/index.php?title=File:GroupAttributesLayout.png> *License:* unknown *Contributors:* ScooterMorris

Image:ForceDirectedSettings.png *Source:* <http://opentutorials.cgl.ucsf.edu/index.php?title=File:ForceDirectedSettings.png> *License:* unknown *Contributors:* ScooterMorris

Image:ForceDirectedDefaultLayout.png *Source:* <http://opentutorials.cgl.ucsf.edu/index.php?title=File:ForceDirectedDefaultLayout.png> *License:* unknown *Contributors:* ScooterMorris

Image:ForceDirected50Layout.png *Source:* <http://opentutorials.cgl.ucsf.edu/index.php?title=File:ForceDirected50Layout.png> *License:* unknown *Contributors:* ScooterMorris

Image:ForceDirected1e-3Layout.png *Source:* <http://opentutorials.cgl.ucsf.edu/index.php?title=File:ForceDirected1e-3Layout.png> *License:* unknown *Contributors:* ScooterMorris

License

Attribution-Noncommercial-Share Alike 3.0 Unported
<http://creativecommons.org/licenses/by-nc-sa/3.0/>